

A composite background image featuring a snowy mountain range, a city skyline, wind turbines, an airplane, and a satellite in the sky.

# ANALYSIS OF SMALL MP AND NP IN SIMPLE AND COMPLEX MATRICES

Andy M. Booth, Lisbet Sørensen

# The 'new' projects



To develop a toolbox of coherent, complementary, analytical instrumental methods and workflows to consistently and reliably detect, identify, quantify (mass/number) and characterize sMP and NP in environmental samples.

**WP3**



Develop a toolbox of complementary analytical methods to reliably detect, identify, and quantify sMP and NP in aquatic environmental samples and to understand their behaviour in different aquatic matrices.

**WP1**




The Research Council  
of Norway






# Microplastic ≠ Nanoparticles

A **200  $\mu\text{m}$**  particle is **10,000x** larger than a **20 nm** particle




# Current status




**Sampling, extraction and analysis approaches for microplastic (>10  $\mu$ m) are not suitable for particles <10  $\mu$ m**

# Challenge #1

How do we go from this?



To this?



Start with simple matrices first!!!

**Plastic size** - small particles are challenging to recover and detect

**Matrix complexity** - high biomass and heterogeneity complicates plastic separation



## BIOLOGICAL SAMPLES

Fish Mussels Plankton




Lack of terrestrial organisms

## PRE-SEPARATION

### DISSECTION

Gastro-intestinal tract or whole body



Recovers >500 μm

## WASTEWATER SAMPLES



Influent


Effluent

Sludge

Lack of proper organic matter content characterization

### IN SITU SIEVING

Plankton net  
Stack filters device  
Mesh size  
10 - 300 μm



Clogging issues  
Loss of smallest fraction

## WATER SAMPLES

Drinking water  
Sea water  
Fresh water



## SEDIMENT SAMPLES

Beach sediments  
Deep sea sediments



## CLEAN-UP & SEPARATION

### DENSITY SEPARATION

#### Flotation

Not suitable for small particles  
Organic matter? Bubbles?

#### Centrifugation


Lack of method development

#### Ultrasonic Separation



### DIGESTION


Plastic degradation?  
pH, temperature?  
Standard protocol?  
Organic matter content?  
Enzymes, H<sub>2</sub>O<sub>2</sub>, HNO<sub>3</sub>, KOH, NaOH, HCl?



### FILTRATION

Filter size  
0.5 - 10 μm

Underestimation of smallest fraction  
because of retention on filter



## ANALYSIS METHODS

### VISUAL

Shape, size, color  
N/A below 200 nm  
Cannot identify chemical composition

#### Microscope SEM, TEM

Plastic spatial structure  
in tissues by microsection

Couple  
with  
staining?



### VIBRATION SPECTROSCOPY

Raman microscopy  
<20 μm

FTIR  
>20 μm or films

↑ Spatial resolution




### MASS SPECTROSCOPY

Suitable for nanoplastics

TDS-GC-MS  
qualitative analysis

py-GC-MS  
greater sensitivity

↑ Sensitivity



**SPIKE CONTROLS** — accurate recovery? effect of plastic size, shape, composition on recovery?

**BLANK CONTROLS** — background level and contamination? False positives from natural organic particles or dyes?

# Where to start?

**To develop and optimise the sample preparation, extraction and concentration steps,  
we first need to establish need robust analysis methods**

**BUT.....**

**.....we need suitable reference materials to develop the analysis methods.**

# Challenge #2

How do we produce a mixture of irregular shaped plastic particles with a size distribution covering 1 nm to 1  $\mu\text{m}$ ?




- Final product can be dry sieved to produce <math><100 \mu\text{m}</math> fraction
- Needs further fractionation using alternative methods

**Cryomilling**



# Reference material goal



# Sample preparation – Filtration

'Standard' filtration approaches do not work for sub-micron particles



**Stirred Cell**



**Cross flow ultrafiltration**

## Challenges

- Filters/membranes that can withstand high pressures.
- Recovery of NMP from filter media.
- Custom-made system?

# Challenge #3

**Which analysis method (s) to use?**



# Quantification – Microplastic (>10 $\mu\text{m}$ )

## Measurements

Polymer identification (individual particles)

Particle counting

Particle size



## Quantification

Total number of particles

Particles of particular polymer type

Particles of particular size (range)




**FPA- $\mu$ FTIR**



**$\mu$ Raman**

# Robust techniques for sMP 1-10 $\mu\text{m}$

- Implementation of automated particle finder in the Raman workflow
- Define instrumental LoDs (size) for  $\mu\text{FTIR}$  and  $\mu\text{Raman}$  as 'start point' for other techniques:
  - Laser confocal Raman - apply autofocus procedures, optimize the laser power and automate workflow.
  - Time-of-Flight secondary ion mass spectrometry (TOF-SIMS) - ion imaging of particles 10-1  $\mu\text{m}$ .
- Investigate combining light or hyperspectral microscopy with Raman spectroscopy and nano-TOF-SIMS.



# Quantification – Nanoplastic (<1 μm)

## Measurements

Polymer identification (individual particles)

Particle counting

Particle size



## Quantification

Total number of plastic particles

Total mass of particular polymer type  
(defined size range)



No single instrument currently capable of achieving this for MP <1 μm and nanoplastic



Possible to achieve



# Techniques for particles $<1 \mu\text{m}$

## Direct Approach

- Optimize cascade filtration down to  $0.45 \mu\text{m}/0.2 \mu\text{m}$  for metal-coated polycarbonate membranes.
- HR SEM as particle finding step.
- Correlative SEM-Raman and SEM-TOF/SIMS to identify polymer particles.



High-Res SEM

# Techniques for particles $<1 \mu\text{m}$



**Py-GC-MS**

## Indirect Approach

Isolate particle fractions  $<0.45 \mu\text{m}$  using cross flow ultrafiltration

Field flow fractionation techniques (Centrifugal or Flow) will be studied

Establish py-GC-MS for plastic quantification in each fraction

# Complex samples – Challenge #4

## SEPARATION METHODS

### DENSITY SEPARATION

#### Flotation

Not suitable for small particles  
Organic matter? Bubbles?

#### Centrifugation


Lack of method development

#### Ultrasonic Separation




### DIGESTION

Plastic degradation?  
pH, temperature?  
Standard protocol?  
Organic matter content?  
Enzymes, H<sub>2</sub>O<sub>2</sub>, HNO<sub>3</sub>, KOH, NaOH, HCl?



### FILTRATION

Filter size  
0.5 - 10 μm




Underestimation of smallest fraction  
because of retention on filter




If the filtration and analysis of sMP and NP is successful for simple matrices, digestion methods for MP can be assessed for sMP and NP.



# Approach we're aiming for



# Nanoplastic in Microred



# Thanks for your attention!

